## LABORATORY EQUIPMENTS







Multiple light chamber for colour comparison



**Standards:** EN 50081-2; EN 60204-1; ASTM D1729; ISO 3668

**ZIPOR** - Equipamentos e Tecnologia Industrial, S.A.

**Commercial Services:** 

**Tel.:** +351.256 831 411 **Fax:** +351.256 831 412

email: pegasil@zipor.com

Technical Services:

**Tel.:** +351.256 888 240 **Fax:** +351.256 888 121

## LABORATORY EQUIPMENTS

Multiple light chamber for colour comparison

## **DESCRIPTION:**

The purpose of the present equipment is to provide controlled conditions for visual analysis in order to ascertain and compare colours.

In structural and generic terms, the present equipment consists of a cabinet of adequate dimensions whose interior has been designed using highly durable materials, and which is matt grey in colour, since this is the way to achieve the most favourable conditions for analysis.

The electrical components of the device, as well as the 5 standard bulbs (**Daylight, Shop Light, Cold White, Incandescent, Ultra-violet**) are installed in the upper part of the cabinet and are duly distributed, always taking into account the safety of the operator. The front command panel houses: 5 push-buttons for activating / deactivating the corresponding bulbs; 1 "POWER" push-button and also an hour counter (with no "Reset" facility) for controlling the operating time vis-à-vis the bulbs, which should be replaced after 2000 hours of use.

Supplied with:

. Operating instructions ; conformity / calibration certificate

Standards:

EN 50081-2; EN 60204-1; ASTM D1729; ISO 3668

**Power Consumption:** 

200 W

Noise:

---- db

**Power Supply:** 

220 V.AC - 50 / 60 Hz

**ZIPOR** - Equipamentos e Tecnologia Industrial, S.A.

Commercial Services: Tel.: +351.256 831 411 Fax: +351.256 831 412

email: pegasil@zipor.com

**Technical Services:** 

**Tel.:** +351.256 888 240 **Fax:** +351.256 888 121